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Адрес: No. 38 Huanzhen Xi Road Dagang Town, Nansha, 511470 Guangzhou, Guangdong,
PEOPLE'S REPUBLIC OF CHINA

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THE MANUFACTURER KEEPS THEIR RIGHTS TO CHANGE SPECIFICATIONS OF THE PRODUCT WITHOUT PRE-WARNING!

IMPORTANT! READ CAREFULLY AND KEEP FOR FUTURE REFERENCE! USE THIS DEVICE EXACTLY AS DESCRIBED IN THE MANUAL AND ONLY IF IT IS IN PROPER WORKING CONDITION.

I. INTENDED USE

The device is an infrared thermometer intended to measure forehead temperature of infants and adults without contact with the human body. It can be used by consumers in household environment and doctor in clinic as reference.

II. RULES AND RECOMMENDATIONS FOR SAFETY USAGE

1. Only use this thermometer once you have read and understood the instruction for use.
2. Only use this thermometer for the specific purpose described in the instruction manual. Any misuse will void the warranty.
3. Children must not play with the device.
4. Device can be used for professional purposes or for personal needs.
5. Using this thermometer is no substitute for consulting a doctor. Do not self-diagnose or self-medicate based on the measurements without consulting your doctor. In particular, do not start taking any new medication or change the type and/ or dosage of any existing medication without prior approval.
6. This thermometer must be operated by adults and be placed out of the reach of infants.
7. **Warning!** Small parts. Risk of swallowing or choking. Make sure that children do not use the device unsupervised.
8. If this thermometer is stored in a cooler or warmer location than where it is being used, put it in the operation place for at least 30 minutes before taking the measurement.
9. After physical activity, a bath, or spending time outdoors, wait approximately 30 minutes before you measure the temperature on your forehead.
10. To ensure an accurate measurement, you should only use this thermometer indoors when the ambient temperature range is between 10°C and 40°C.
11. Before and after each use, the thermometer should be cleaned according to the requirement of this manual to avoid possible bacterial infections.
12. Do not expose this thermometer to electric shocks.
13. Do not expose the thermometer on extreme temperature conditions, bigger than 55° C and less than - 20° C.
14. Do not use the device in relative humidity, higher than 85%.
15. Do not put this product in a place with direct sunlight, high temperature, high humidity, excessive dust, or the place where it may contact corrosive gas or liquid. Do not use this product in the above environment.
16. Never drop the device.
17. **DO NOT TOUCH THE GLASS OF THE INFRARED SENSOR WITH YOUR FINGERS AS IT IS EXTREMELY FRAGILE AND COULD CAUSE INJURY.**
18. Do not use this thermometer if it is damaged.
19. Do not open up this thermometer. It contains small parts which might be swallowed by children.
20. Never immerse the device in water or other cleaning liquids, and do not let any liquid enter the interior of this product.
21. If you will not use the device for long time, take the batteries out from the battery compartment. This way you will prevent the device from battery leaking or corrosion.
22. Keep the thermometer away from electromagnetic fields produced by objects such as microwave ovens, arc welders or induction heaters.
23. The device should no longer be reused when it reaches the expected service life (5 years, from the date of production) of the device. See chapter XIII to dispose the device.
24. Electrical components and batteries have heavy metal or chemicals which may do harm to environment, follow the local law and regulation on its disposal.
25. If a problem with the device appears, contact with the trader, from which you have bought it. Do not try to repair the device by yourself. Otherwise the warranty will be void:
 - If an error occurs or if the thermometer is damaged, DO NOT attempt to repair the thermometer yourself, as this will void the warranty. Only have repairs carried out by authorized service partners.
 - If you need assistance in setting up, using, or maintaining the thermometer.
 - If you have unexpected operation or events to report.

III. WORKING PRINCIPLE

All objects, solid, liquid or gas emit energy by radiation. The intensity of this energy depends of the temperature in the object. The infrared thermometer Easy Check is therefore able to measure the human's temperature, only by the energy that person emits. This measurement is possible, thanks to an external temperature probe on the device, which permanently analyses and register the temperature of the environment.

Therefore, immediately after the thermometer is closed to the body and the radiation sensor is activated, the temperature is immediately measured through detecting infrared heat, caused by arterial blood flow. The temperature of the body can be measured also without any interference from the heating of the surrounding environment.

Note: The reference body temperature for the non-contact infrared thermometer is axillary temperature.

IV. DEVICE DESCRIPTION

➤ **See Figure PD:**

1. Display; 2. MODE button; 3. Button to change the unit (Celsius or Fahrenheit); 4. Memory button; 5. Infrared sensor; 6. ON/ SCAN button; 7. Cover of the battery compartment.

➤ **See Figure D:**

8. Celsius; 9. Measured temperature; 10. Fahrenheit; 11. Memory symbol; 12. Low battery symbol; 13. Buzzer/ alarm symbol.

V. INSTRUCTIONS FOR USING

1. Put the batteries in the compartment correctly according its polarity. See Figure B.

2. If using the device for the first time or when adding new batteries, wait for 10 to 15 minutes for the device to warm up. This will allow to the device to acclimate to the room temperature.

3. Press the on/scan button, point the device to the forehead (look diagram for pointing the thermometer), on distance from 3 to 5 cm. Press the button ON/ Scan in standby mode and the measuring will finish when the temperature shows itself on the display or a sound is issued. Measuring time is one second. **Note: Do not change the position of the thermometer before measuring to be finished.**

4. Before temperature measuring, ensure that you have removed hair or sweat from the forehead.

➤ **How to measure temperature**

Hold the thermometer against forehead, in the right area, on distance from 3 to 5 cm. Press the Scan button for measuring the temperature and it will immediately appear on the display.

For reliable results follow the guidelines below:

1. Do NOT attempt to take an individual's temperature by aiming the device at a body part other than the forehead (i.e., arms, legs, torso, etc.).

2. Stay still and do NOT move the thermometer before the final result is displayed.

3. Do NOT take an individual's temperature near any drafts (i.e. wind, air conditioner, fans).

4. Do NOT take an individual's temperature shortly after he or she participated in vigorous physical activity, bathed, took medicine, ate or drank, or had a cold or warm compress on his or her forehead.

5. Do NOT take an individual's temperature when there has been a significant change in ambient conditions. Allow the device to adjust to the significant change in the ambient temperature due to a change in environment for at least 30 minutes before taking the measurement.

6. Wait at least 3 to 5 seconds between measurements.

7. Remove any hair from the forehead.

8. Wipe the forehead, if there is any sweat on it.

9. We recommend that you take three measurements on the forehead. If the results of the three measurements are different, use the one with the highest value.

VI. MENU SETTINGS AND FUNCTIONS

1. Turning the device on:

Press the ON/ Scan button, after a second, on the display will be shown symbols „—°C“ or „—°F“. Then press the ON/ SCAN button once more and the result will show up automatically after 1 second. If you do not use the thermometer anymore, it will turn off automatically after 30 seconds.

2. Mode switching

a. Press the button "MODE" and the display will show: Body°C

b. Press the button "MODE" again and display will show : Room ...°C

c. Press again the button "MODE" and display will show: Surface Temp°C

Note: The thermometer is set by default on „BODY“ mode.

IMPORTANT!

Surface temperature differs from the internal temperature of human body. For measuring the body temperature, always use „BODY“ mode. Make sure you select „SURFACE TEMP“ mode to read the external area only.

3. In turned on state, press the button, "C/F", to switch between Celsius and Fahrenheit.

4. In turned on state, press the button "MEM" (Memory), which will show the last measured temperature and will show last 32 measurements. If again, in turned on state, you press the button "MEM" and hold it for 5 seconds, all data will be deleted. After this press again „MEM“ button and the display will show „CLr“.

5. In turned on state, press the buttons "MODE" and "C/F" at the same time for 3 seconds. The display will show: "F4", and then press "MEM" button to enter into F5 (ON or OFF functions). Press MODE to display "ON" and the buzzer is active. Press "C/F" to display "OFF", and the buzzer will not be active. Press "MEM" to confirm setting.

6. Device calibration through menu F4: When there is difference between Easy Check and mercury thermometer, you believe that the mercury thermometer is more accurate but not comfortable for using. You can use calibration function to tune Easy Check and to receive the value, same as this on the mercury thermometer. Furthermore, when you use Easy Check for people with different skin color (for example from yellow, white or black race), you can also use calibration.

> Calibration instructions

a. In turned on state, press the buttons "MODE" and "C/F" at the same time for 3 seconds.

b. The display will show: "F4".

c. Press "MODE", to increase with 0.1°C, press "C/F", to decrease with 0.1°C.

d. Press "MEM" button twice, to keep the settings.

e. In case of seasonal or ecological changes, a check and device regulation must be done.

VII. TECHNICAL CHARACTERISTICS

Items	Specifications
Battery	DC 3V (2 pcs AAA batteries)
Temperature sensor	Infrared Sensor
Measuring distance	3cm - 5cm (1.2in - 2in)
Measurement range	Body temperature mode: 32.0°C ~ 43.0°C (89.6°F ~ 109.4°F) Under body mode, there is three color backlights: Green color backlit: ≤ 37.3°C (99.1°F), means normal temperature. Orange color backlit: 37.4°C - 37.9°C (99.3°F - 100.2°F), means may be low fever. Red color backlit: ≥ 38°C (100.4°F), means may be high fever. Surface temperature mode: 0°C ~ 60°C (32°F ~ 140°F) Room temperature mode: 0°C ~ 40°C (32°F ~ 104°F)
Measuring accuracy	Body temperature mode: 32.0°C ~ 34.9°C (89.6°F ~ 94.8°F) : ±0.3°C (±0.6°F) 35.0°C ~ 42.0°C (95°F ~ 107.6°F) : ±0.2°C (±0.4°F) 42.1°C ~ 43.0°C (107.8°F ~ 109.4°F) : ±0.3°C (±0.6°F)
Clinical repeatability	≤ 0.2°C (0.4°F)
Display resolution	0.1°C (0.1°F)
Measurement Time	1 second
Memory	32 sets
Operation Environment	Ambient temperature: 10°C ~ 40°C (50°F ~ 104°F) Relative humidity: ≤85%RH Atmospheric Pressure: 700 hPa to 1060 hPa
Storage and transportation condition	Ambient temperature: -20°C ~ 55°C (-4°F ~ 131°F) Relative humidity: ≤ 95%RH Atmospheric Pressure: 700 hPa to 1060 hPa
Weight	Approx. 67g (Battery excluded)
Product Size	Approx. 130 x 45 x 55 mm (L x W x H)
Expected Service Life	5 years
Automatic power-off	≤ 30 secs

Items	Specifications
Measuring site	Forehead
Reference body site	Axillary
Mode of operation	Adjusted mode

VIII. BATTERIES

1. Battery change:

a. When LCD display shows this symbol , it means the batteries are exhausted and should be changed with new ones.
b. To change the batteries with new ones, open the cover and remove the exhausted batteries and place the new ones (See Figure B). Be careful to follow the correct polarity when placing the new batteries. An error with this action can lead to damage of the device and void of the warranty. Never use rechargeable batteries.

2. DISPOSING OF BATTERIES: To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations. The disposal method is defined in Chapter XIII.

IMPORTANT BATTERY GUIDELINES

➤ IMPORTANT BATTERY GUIDELINES:

1. An adult must handle the battery.
2. Do not allow children to play with the batteries. Risk of swallowing! Small children can swallow the batteries and to suffocate. Store the batteries at place where children cannot reach!
3. Always remove exhausted batteries and do not throw them away together with the household waste, but in the intended places. They are recyclable.
4. Use batteries only AAA size.
5. Alkaline batteries are recommended.
6. When you change with new batteries, always change all batteries. Children must not be present when changing the batteries.
7. Always place the batteries correctly according to their polarity (+/-), marked on the bottom of the battery compartment. If they are not correctly placed, this can lead to short circuit or electric shock.
8. Do not mix old batteries with new ones.
9. Do not mix alkaline, standard (Carbon-Zinc) or rechargeable (Nickel- Cadmium) batteries.
10. Do not put the batteries in fire, because they can explode or to leak.
11. Keep the batteries from direct heating sources.
12. If you will not use the device for a long time, take the batteries out of the battery compartment. This way you will prevent the device from battery leaking or corrosion.
13. If when changing the batteries, you find that the batteries leak, put safety gloves on and clean the battery compartment with dry towel.
14. Non-rechargeable batteries must not be charged.
15. Do not disassemble, open or shatter the batteries.
16. If your skin or eyes come into contact with battery fluid, rinse the affected area with water and seek medical assistance.
17. Protect batteries from excessive heat.

IX. INCLUDED IN THE PACKAGE

- Non-contact Infrared Thermometer - 1 pc.
- User Manual - 1 pc.
- AAA alkaline batteries - 2 pcs.
- Storage case - 1 pc.

X. CLEANING AND PRODUCT MAINTENANCE

1. The thermometer is a multi-patient reusable device, and the cleaning and disinfection instructions were validated. Follow the instructions to clean/disinfect your thermometer prior to use.

2. Cleaning:

- a. Use a cotton swab moistened with hot water or neutral detergent to throughout clean the casing of thermometer and measuring probe to remove stain or oil.
- b. Inspect the thermometer visually after the cleaning process, ensure there is no soil or stain is seen. If the device is determined not to be visually clean at the end of the cleaning step, you should either repeat the relevant previous cleaning steps or safely dispose of the device so that a visibly soiled device is not used again.

3. Disinfection:

- a. Disinfect the thermometer casing and the probe for 30 seconds with a cotton swab or non-dust cloth dipped in 75% ethanol.
- b. Repeat the above steps at least 3 times, then wipe surface and probe dry with a clean cloth, and wait about 1 minute for drying after disinfection process.

4. Note:

- Ensure that no liquid enters the interior of the thermometer.
- Never use abrasive cleaning agents, thinners, or benzene for cleaning, and never immerse the thermometer in water or other cleaning liquids.
- Take care not to scratch the surface of the probe lens and the display.
- If the visual inspection identifies unacceptable deterioration such as corrosion, discoloration, pitting, cracked, the device should not be further used and should be discarded.
- If the thermometer has been used in patients, visually inspect the thermometer to determine if there is stain or oil on the casing or probe after each use, and follow the instruction of cleaning and disinfection to cleaning your thermometer.

XI. TROUBLESHOOTING

If you have problems, while using the thermometer, take a look at the guidelines in this manual for solving the problem. If the problem still exists, contact with the trade agent, which you bought the product from.

1. Display shows temperature, higher than 43°C (109.4°F)

The temperature is in Fahrenheit. Change the measurement to Celsius.

2. Display shows temperature lower than 32 °C (89,6 °F)

For measuring the surface temperature, press button "BODY" and set for measurement on „BODY“. If the device is in mode „Surface temperature“, the temperature 32 °C (89,6 °F), which display will show, taking into account the external body temperature, instead internal.

3. Display show the message "Hi"

When using contactless thermometer, model Easy Check, the message "Hi" can show on the display. In this case the temperature is either over the chosen range or over 43°C (109.4°F) in „BODY“ mode.

4. Display shows message "Lo"

When using the contactless thermometer, model Easy Check, the message "Lo" can show on the display. In this case, recorded temperature is either under chosen range or under 32 °C (89,6 °F) in „BODY“ mode.

The display shows this message because of many reasons. Check the below list with the main possible problems:

Reasons the display shows the message "Lo"	Advise
Temperature measurement is difficult from hair or sweating.	Ensure that there are no obstacles or moisture before temperature measurement.
Temperature measurement is difficult from air flows or sharp changes in environment temperature.	Ensure that there are no air flows around the area, in which the temperature is measured, this can affect the infrared reading.
Successive measurements of the temperature are too close one after another and the thermometer do not have time to restart.	Leave interval 3-5 secs. Between measurements; recommended is pause for 15 sec.
Measurement area is away from the thermometer.	Measure the temperature on recommended distance (around 3-5 cm, 1.2 ~2.0 inches).

XII. SYMBOL EXPLANATION

Symbol	Explanation
	Trade mark
	Indicates the need for the user to consult the instructions for use for important cautionary information such as warning and precautions that can not, for a variety of reasons, be presented on the medical device itself.
	IEC 60417-5333, type attached part BF
	IEC 60417-5031, Direct current

Symbol	Explanation
IP22	Degree of protection against harmful ingress of water and particulate matter 2 Protected and against solid foreign objects of 12.5mm Φ and greater 2 Protection against vertically falling water drops when ENCLOSURE tilted up to 15°
	Refer to instruction manual
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Specific serial number
	Lot number
	Importer
	Model (Item) number
	Medical device
	Manufacturer
	Country of manufacture
	Symbol of European Authorized Representative, expressly indicate the name and address of European Authorized Representative after the symbol.
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical device.
	Date of Manufacture
	Unique device identifier

XIII. DISPOSAL

It is required that you should not dispose of WEEE as unsorted municipal waste and collect such WEEE separately. In accordance with Directive 2012/19/EU and national disposal regulations regarding old electrical and electronic

devices, be advised that such items must be disposed of in a special way within the European Union (EU). These regulations require the environmentally friendly recycling/disposal of old electrical and electronic devices. Such items must not be disposed of as domestic refuse. This has been expressed using the icon of the "crossed out trash can".

The Non-contact Infrared Thermometer (Model: JXB-183) contains batteries and recyclable electronic waste. To protect the environment, do not dispose of it in household waste as unsorted municipal waste, but collect separately and take it to an appropriate local collection center. The Non-contact Infrared Thermometer (Model: JXB-183) is classified as small equipment (5th category) according to ANNEX IV of Directive 2012/19/EU.

Before disposing of your Non-contact Infrared Thermometer (Model: JXB-183), take away the batteries from Non-contact Infrared Thermometer (Model: JXB-183) and dispose of the batteries and Non-contact Infrared Thermometer (Model: JXB-183) separately.

You can take away the batteries as follows – See Figure C:

- a. Open the lid on the back of the device.
- b. Remove the batteries.

Batteries represent a hazard to the health and the environment! Never open, damage, or swallow batteries or allow them to pollute the environment. They may contain toxic, ecologically hazardous heavy metals. The removed batteries shall be disposed of at the point of sale or in the corresponding containers provided at collection points by local public waste authorities and treated according to 2006/66/EC.

The batteries and recyclable electronic waste shall be disposed of in a separate waste container. Packaging materials must be disposed of according to local regulations. Consult the authorized collection points for more information.

In Germany: For the recyclable electronic waste disposal, you can search the authorized collection points at the following website: <https://www.ear-system.de/ear-verzeichnis/sammel-und-ruecknahmestellen>. For batteries, you can search the authorized collection points at the following website: <https://www.ear-system.de/ear-verzeichnis/battghersteller#no-back>. Search method: using your postcode or the name of the city and state where you live to find collection points near you.

Other countries: For country-specific information on disposal, contact your local dealers or importer.

XIV. DECLARATION OF CONFORMITY

The appliance is suitable for home health care and clinic environments.

Warning! Do not use near high frequency surgical equipment or in a RF shielded room of an ME system for magnetic resonance imaging, where the intensity of the EM disturbances is high.

Warning! The use of this equipment near to or together with another equipment should be avoided, because this may result in improper operation. If such use is necessary, both of the equipments should be observed to verify that they are operating normally.

Warning! The use of accessories, transducers and cables other than those specified or provided by the manufacturer may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning! Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment JXB-183, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could occur.

Follow all instructions to ensure the basic safety and essential performance of the appliance. The portable or mobile RF communication equipment may affect the performance of the JXB-183. Avoid using the appliance near mobile phones, microwave ovens or other appliance, which emit strong electromagnetic waves.

GUIDELINES AND MANUFACTURER'S DECLARATION OF THE ELECTROMAGNETIC EMISSIONS AND IMMUNITY

TABLE 1:

Guidance and Manufacturer's Declaration – Electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	GROUP 1
RF emissions CISPR 11	CLASS B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage Fluctuations/ Flicker emissions IEC 61000-3-3	Not applicable

TABLE 2:

Guidance and Manufacturer's Declaration – Electromagnetic emissions		
Immunity test	IEC 606011-2 level test	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 6 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/ burst IEC 61000-4-4	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable

Guidance and Manufacturer's Declaration – Electromagnetic emissions		
Immunity test	IEC 606011-2 level test	Compliance level
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable
Power Frequency Magnetic field IEC 61000-4-8	30 A/m 50Hz/ 60Hz	30 A/m 50Hz/ 60Hz
Conducted RF IEC 61000-4-6	Not applicable	Not applicable
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz
Note: U_T is the a.c. mains voltage prior to application of the test level.		

TABLE 3:

Guidance and Manufacturer's Declaration – Electromagnetic immunity						
Radiated RF IEC 61000-4-3 (test specifications for Enclosure Port Immunity to RF wireless communications' equipment)						
Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation	Distance (m)	Immunity test level (V/M)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710 745 780	704-787	LTE Band 13,17	Pulse modulation 217 Hz	0.2	0.3	9
810 870 930	800-960	GSM800/900, TETRA800, CDMA850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
1720 1845 1970	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25;UMTS	Pulse modulation 217 Hz	2	0.3	28
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240 5500 5785	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9

This appliance is manufactured by: Guangzhou Berrcom Medical Device Co., Ltd
Address: No. 38 Huanzhen Xi Road Dagang Town, Nansha, 511470 Guangzhou, Guangdong,
PEOPLE'S REPUBLIC OF CHINA

The appliance is distributed in Bulgaria by: Moni Trade Ltd; **Address:** Bulgaria, Sofia, Trebich, Dolo str. 1 **Tel.:** +359 2/ 936 07 90
Authorized representative of the manufacturer for Europe: Shanghai International Holding Corp. GmbH (Europe)
Address: Eiffestrasse 80, 20537 Hamburg, Germany

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Acest aparat este fabricat de: Guangzhou Berrcom Medical Device Co., Ltd.
Aparatul este distribuit în Bulgaria de: Moni Trade Ltd; **Adresa:** Bulgaria, Sofia, Trebich, Dolo str. 1