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Section 1: General Information Section 2: Project Details Section 3: Financials Section 4: Risk Assessment Section 5: Compliance Section 6: Other Information									
Project Name		Project ID		Project Manager		Project Status		Project Budget	
Project A		P001		John Doe		On Track		\$100,000	
Project B		P002		Jane Smith		Delayed		\$200,000	
Project C		P003		Mike Johnson		On Hold		\$50,000	
Project D		P004		Sarah Brown		Completed		\$75,000	
Project E		P005		David White		On Track		\$150,000	
Project F		P006		Emily Green		On Track		\$120,000	
Project G		P007		Chris Black		On Track		\$90,000	
Project H		P008		Alex Blue		On Track		\$110,000	
Project I		P009		Mia Purple		On Track		\$130,000	
Project J		P010		Noah Yellow		On Track		\$140,000	

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Shenzhen AOJ Medical Technology Co., Ltd.

Add: Room 301&4F, Block A, Building A, Jingfa Intelligent Manufacturing Park, Xiawei-yuan, Gushu Community, Xixiang Street, Bao'an District, 518126 Shenzhen, CHINA

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THE MANUFACTURER KEEPS THEIR RIGHTS TO CHANGE SPECIFICATIONS OF THE PRODUCT WITHOUT PRE-WARNING!

IMPORTANT! READ CAREFULLY AND KEEP FOR FUTURE REFERENCE! USE THIS DEVICE EXACTLY AS DESCRIBED IN THE MANUAL AND ONLY IF IT IS IN PROPER WORKING CONDITION.

I. INTENDED USE

The appliance is infrared thermometer is intended for measuring the body temperature of both children and adults. The temperature can be measured either by pointing the thermometer at the child's or adult's forehead or by auricular measurement (outer ear canal). Measuring the temperature by the auricular method is suitable for children above 3 months of age. The device can be used in home care environment. The product is tested according to and complies with the requirements of ISO 81060-2-56, ISO10993-5, ISO10993-10 and ISO10993-23.

II. INCLUDED IN THE PACKAGE

Before use, carefully open the package and check that all parts are present as listed below and have not been damaged during shipment. If there is any damage or problem related to operation, contact the dealer from whom you purchased the product.

INCLUDED IN THE PACKAGE:

Infrared thermometer – 1 pc;

Pouch – 1 pc;

User manual – 1 pc;

Batteries AAA 1.5V – 2 pcs.

III. RULES AND RECOMMENDATIONS FOR SAFE USAGE

1. Only use this thermometer once you have read and understood the instruction for use.
2. Only use this thermometer for the specific purpose described in the instruction manual. Any misuse will void the warranty.
3. Do not place the thermometer near the reach of children as this may lead to damaging the product.
4. Children must not play with the device.
5. Device can be used for personal needs.
6. Do not use the thermometer to measure the temperature of newborn babies or for continuous temperature monitoring purposes.
7. Using this thermometer is no substitute for consulting a doctor. Do not self-diagnose or self-medicate based on the measurements without consulting your doctor. In particular, do not start taking any new medication or change the type and/ or dosage of any existing medication without prior approval. Consult your doctor if you have symptoms such as vomiting, diarrhea, dehydration, loss of appetite and fatigue, seizures, shivering, stiff neck, pain when urinating, etc., even in the absence of fever.
8. This thermometer must be operated by adults and be placed out of the reach of infants.
9. Warning! Small parts. Risk of swallowing or choking. Make sure that children do not use the device unsupervised.

- 10.** If this thermometer is stored in a cooler or warmer location than where it is being used, put it in the operation place for at least 30 minutes before taking the measurement.
- 11.** After physical activity, a bath, or spending time outdoors, wait approximately 30 minutes before you measure the temperature on your forehead.
- 12.** To ensure an accurate measurement, you should only use this thermometer indoors when the ambient temperature range is between 10°C and 40°C
- 13.** Before and after each use, the thermometer should be cleaned according to the requirement of this manual to avoid possible bacterial infections.
- 14.** Do not forcibly insert the probe into the ear canal, otherwise you may injure the ear canal. If the ear canal is dirty, first clean it with cotton swab and then measure the temperature.
- 15.** Do not measure the - Do not measure your temperature auricularly if you have an ear infection - otitis.
- 16.** Do not expose this thermometer to electric shocks.
- 17.** Do not expose the thermometer on extreme temperature conditions, bigger than 40°C and less than 10°C.
- 18.** Do not use the device in relative humidity, higher than 95%.
- 19.** Do not put this product in a place with direct sunlight, high temperature, high humidity, excessive dust, or the place where it may contact corrosive gas or liquid. Do not use this product in the above environment.
- 20.** Never drop the device.
- 21.** DO NOT TOUCH THE GLASS OF THE INFRARED SENSOR WITH YOUR FINGERS AS IT IS EXTREMELY FRAGILE AND COULD CAUSE INJURY.
- 22.** Do not use this thermometer if it is damaged.
- 23.** Do not open up this thermometer. It contains small parts which might be swallowed by children.
- 24.** Never immerse the device in water or other cleaning liquids, and do not let any liquid enters the interior of this product.
- 25.** If you will not use the device for a time longer than 2 months, take the batteries out from the battery compartment. This way you will prevent the device from battery leaking or corrosion.
- 26.** Keep the thermometer away from electromagnetic fields produced by objects such as microwave ovens, arc welders or induction heaters.
- 27.** Portable and mobile radio frequency communication devices may affect the operation of the infrared thermometer.
- 28.** The device should no longer be reused when it reaches the expected service life (5 years, from the date of production) of the device. Electrical components and batteries have heavy metal or chemicals which may do harm to environment, follow the local law and regulation on its disposal.
- 29.** If a problem with the device appears, contact with the shop, from which you have bought it. Do not try to repair the device by yourself. Otherwise, the warranty will be void:
- If an error occurs or if the thermometer is damaged. DO NOT attempt to repair the thermometer yourself, as this will void the warranty.
 - Only have repairs carried out by authorized service partners.
 - If you need assistance in setting up, using, or maintaining the thermometer.
 - If you have unexpected operation or events to report.

30. The replaceable parts of the thermometer are the 2 batteries AAA 1.5V. Replace them as soon as they have no charge. Do not use batteries with different voltages or different specifications. The product does not have any other replaceable parts and materials.

IV. WORKING PRINCIPLE

All objects, solid, liquid or gas emit energy by radiation. The intensity of this energy depends of the temperature in the object. The infrared thermometer Easy Check is therefore able to measure the human's temperature, only by the energy that person emits. This measurement is possible, thanks to an external temperature probe on the device, which permanently analyses and register the temperature of the environment. Therefore, immediately after the thermometer is closed to the body and the radiation sensor is activated, the temperature is immediately measured through detecting infrared heat, caused by arterial blood flow. The temperature of the body can be measured also without any interference from the heating of the surrounding environment. The temperature can be measured and by placing the sensor in the ear. The product is with alarm for high temperature, front light, mute button, function to change the measurement unit °C and °F, automatic turning off to save energy.

V. DEVICE DESCRIPTION


FIGURE A: Device

1. LED display screen
2. Mute-unmute button/Unit toggle button (°C, °F) .
3. Measure button/ Power on/off button/Front light on/off button
4. Mode button (human body/adjusted temperature)/Memory
5. LED indicator light

6. Probe cover (put the cover on when measuring forehead temperature)
7. Probe (take the cover off when measuring ear temperature)
8. Front light (it lights up in ear temperature mode)
9. Battery compartment
10. Battery compartment button

FIGURE B: Display


1. Adult forehead temperature mode (above 12 years old)
2. Child forehead temperature mode
3. Ear temperature mode
4. Adjusted mode
5. Fahrenheit degree °F
6. Celsius degree °C
7. Temperature value
8. Front light indicator
9. Low battery indicator
10. Mute/Unmute indicator

VI. Sound and LED indicator description - TABLE C:


Temperature range	Sound	LED Indicator status
Forehead temperature		
32.0°C-37.5°C / 89.6°F - 99.5°F	One long beep	Green
37.6°C - 38.0°C / 99.6°F -100.4°F	6 short beeps (2 short beeps each time)	Orange
38.1°C - 42.9°C / 100.5°F - 109.2°F	6 short beeps (2 short beeps each time)	Red
Ear temperature (external ear canal)		
32.0°C-37.5°C / 89.6°F - 99.5°F	One long beep	Green
37.6°C - 38.0°C / 99.6°F -100.4°F	6 short beeps (2 short beeps each time)	Orange
38.1°C - 42.9°C / 100.5°F - 109.2°F	6 short beeps (2 short beeps each time)	Red
Adjusted temperature (object temperature)		
0°C-100°C / 32.0°F -199.0°F	One long beep	

VII. DISPLAY AND OPERATIONS DESCRIPTION/MENU SETTING – TABLE D


SCREEN DISPLAY	OPERATING INSTRUCTIONS/MENU SETTING	SOUND AND LED INDICATOR STATUS
MEASURING EAR TEMPERATURE		
	Take off the probe (7) cover (6), press the measure button (3) to turn on the thermometer, and the thermometer will automatically enter ear temperature mode and the icon will display on the screen. Or you can also press mode button (M button) (4) to switch to ear temperature mode if necessary.	See Table C
	Properly insert the probe into the ear canal, press the measure button (3), and the ear temperature will display on the screen.	See Table C
MEASURING FOREHEAD TEMPERATURE		
	Put on the probe (7) cover (6), press the measure button (3) to turn on the thermometer, then press the mode button (M button) (4) to switch to human body temperature mode, and the thermometer will enter forehead temperature mode and the icon will display on the screen.	See Table C
	Aim the probe at the center of the forehead (between the eyebrows) until touching the forehead skin or with a distance of 0-3cm (0-1.18 inches), press the measure button (3), and the forehead temperature will display on the screen.	See Table C
DISPLAY SHOWS THE MESSAGE "HI"		
	When using contactless thermometer the message "Hi" can show on the display. In this case the temperature is either over the chosen range or over 42,9 °C (109,2 °F) in „Ear temperature“ or „Forehead temperature“ mode.	One long beep and LED indicator in red

DISPLAY SHOWS MESSAGE "LO"		
	When using the contactless thermometer, the message "Lo" can show on the display. In this case, recorded temperature is either under chosen range or under 32 °C (89,6 °F) in „Ear temperature“ or „Forehead temperature“ mode.	One long beep and LED indicator in red
CHECKING 40 SETS OF MEMORY DATA – MEMORY FUNCTION		
	In power-on status or power-off status, press and hold the memory button (M button) for about 5s to enter memory mode. The screen display "01", then displays the 1st group of recorded temperature value. 01 indicates the last recorded temperature value. Press the memory button again to check the next recorded value, and the screen displays "02", then displays the recorded value. Continue to press the Memory button to view the next previously stored measurement. You can check a total of 40 temperature readings. Any new measurement will be recorded and the oldest memory deleted without you having to do anything.	No sound

SWITCHING BETWEEN MUTE/ UNMUTE

	In power-on status, press the mute/unmute button  (2) to turn on or turn off the sound. When the sound is turned on, there will be a long beep; when the sound is turned off, the mute icon  will be displayed.	One long beep when the sound is turned on.
--	--	---

SWITCHING BETWEEN ADJUSTED/FOREHEAD/EAR TEMPERATURE MODE

	Press the mode button (M button) (4) to switch between adjusted temperature  mode and human body temperature modes.	No sound
	Human body temperature mode includes: forehead temperature ( or ) mode and ear temperature () mode.	No sound

SWITCHING BETWEEN FAHRENHEIT/ CELSIUS

	In power-on status, press and hold the unit toggle button (2) for about 5s to switch between °C and °F. The selected system will appear on the display (°C and °F). By holding the button (2) again for 5 seconds, you will switch to the other system.	No sound
--	--	-----------------

ERROR MESSAGE/ LOW BATTERY		
	When using the contactless thermometer, the message "Er1" can show on the display. In this case, the ambient temperature is higher than 40.0°C (104.0°F) or lower than 0.0°C (32.0°F).	One long beep and LED indicator in red
	When using the contactless thermometer, the message "ErC" can show on the display. In this case, internal storage data error or temperature calibration is not completed. The thermometer will automatically turn off after 3s.	One long beep and LED indicator in red
 	When the thermometer is powered on and activated and the battery voltage is between 2.4V and 2.7V, then low battery icon will flash. However, this will not affect your normal use.	No sound
	When the thermometer is powered on and activated and the battery voltage is lower than 2.4±0.1V, then only the low battery icon will be displayed. The thermometer will automatically turn off after 3s.	No sound
Manual power-off	In non-ear temperature mode, press and hold the measure button (3) for 5s, and the thermometer will turn off.	The screen turns off

Note: If the battery icon starts flickering, this means that the batteries have no charge. Replace them with new ones to avoid incorrect results.

VIII. METHODS OF TEMPERATURE MEASUREMENT

A. MEASURING BODY TEMPERATURE IN THE EAR (FOR AGE OVER 3 MONTHS)

E

1. Before first time use, place batteries in the device.
2. Remove the probe (7) cover (6).
3. Press the Power/Measure button (3) to turn on the device. The thermometer will enter into ear mode automatically (see Table D). You can also press the mode button M(4) and choose ear measuring mode.
4. Place the probe in the ear as show in figure E. Press the Power/Measure button. The measured temperature will be displayed on the screen immediately.
FOR CHILDREN UNDER ONE YEAR OLD: As the ear canal is slightly curved, you have to pull the ear slightly backwards before inserting the sensor tip. This is important so that the sensor tip can be pointed directly at the ear drum.
FOR CHILDREN ABOVE ONE YEAR OLD AND ADULTS: As the ear canal is slightly curved, you have to pull the ear slightly up and backwards before inserting the sensor tip. This is important so that the sensor tip can be pointed directly at the ear drum.

WARNING! Do not forcibly insert the probe into the ear canal; otherwise, you may injure the ear canal.

WARNING! Be careful when measuring the ear temperature of a child, as a child's earhole is small.


B. MEASURING BODY TEMPERATURE IN FOREHEAD AREA

Practical considerations when taking a temperature

WARNING!

- To achieve optimal measurement accuracy, the forehead area should be dry, clean of sweat and cosmetics. Taking vasoconstrictive medications and skin irritations on the forehead can distort the result when measuring forehead temperature.
- Keep in mind that forehead temperature can be affected by prolonged exposure to the outdoors (e.g. cold air in winter or exposure to the sun in summer) as well as by wearing a hat.
- Forehead temperature measurement is also affected by the correct placement of the thermometer. The thermometer should be placed against the center of the forehead at a distance of 0-3 cm (0-1.18 inch) from the skin surface (see Figure F).

FIGURE F: Correct placement of the forehead thermometer

F

1. If you use the thermometer for the first time, make sure that the batteries are placed. The cover (6) should be in position on the probe nozzle (7). Press the Power/Measure button (3) to turn on the device.
2. Press the mode button (4) to choose forehead mode (adult / child): child temperature measurement mode and adult measurement mode (over 12 years of age) (See Table D for more information).
3. Point the probe at the center of the forehead (between the eyebrows) with a distance of 0-3cm (0-1.18 inches).
4. Press the Power/Measure button. The measured temperature will be displayed immediately on the screen.
5. The thermometer will automatically turn off in 30 seconds after you stop using it.

C. ADJUSTED TEMPERATURE MEASURING (ROOM/OBJECT TEMPERATURE)

G

1. The cover (6) should be in position on the probe nozzle (7). Press the Power/Measure button (3) to turn on the device.
2. Touch the mode button (4) to choose adjusted mode.
3. Point the temperature probe at the adjusted objects with a distance of 0-3cm(0-1.18 inches)
4. Press the Power/Measure button. The measured temperature will be displayed on the screen. The thermometer will automatically turn off in 30 seconds after you stop using it.

D. After measuring the temperature:

1. After each measurement you can switch to memory mode to memorize the measured temperature. See table D.
2. After you stop using the product, clean it and store it in a dry and well-ventilated place.
3. You should wait at least 10 seconds between each measurement.

Note:

1. The thermometer is intended for indoor use only. Ensure that during using the thermometer there is no strong air convection (for example, from electric fan, air

conditioner or a heater).

2. Ensure that the ear canal is clean before measuring the ear temperature. If there is earwax, it is recommended to clean it with cotton swab to avoid incorrect results and contamination of the sensor.

3. Do not hold the thermometer for too long as it is sensitive to the ambient temperature.

4. Make sure that the sensor is clean before use.

5. When measuring the forehead temperature, make sure that the forehead is dry and not sweaty, and there is no hair above it. If the forehead is sweaty or before measuring the temperature you have applied cold compresses, this will lead to incorrect results.

IX. PRODUCT TECHNICAL SPECIFICATIONS

Model:	AOJ-20F
Battery-powered mode/ Batteries	DC 3V/ 2 batteries 1.5 V AAA
Display mode:	LED
Display temperature resolution:	0.1°C/0.1°F
Working mode:	Continuous operation
Measurement time:	Short recall time - around 1 second
Measuring range:	Forehead/ Ear temperature: 32°C – 42.9°C (89.6°F – 109.2°F) ±0,2°C / ±0,4°F.
	Adjusted temperature (room/object temperature): 0.0°C-100.0°C (32.0°F-199.0°F) ±1.0°C/±2.0°F.
Display resolution:	0.1°C
Automatic shutdown:	The thermometer will shut down automatically after 30 seconds if is not being used.
Memory data:	40 measurements
Product size:	16.5 cm x 4.0 cm x 4.3 cm
Weight (g):	about 75 g (without batteries)
Normal operating conditions (Using environment):	- Environment temperature: from +10°C to +40°C (+50°F - +104°F) - Relative humidity: 15%~95% RH - Atmospheric pressure: 70kPa~106kPa
Storage and transportation condition (storage environment):	- Environment temperature: -20°C~+55°C (-4°F-131°F) - Relative humidity: 0~95% RH - Atmospheric pressure: 50kPa~106kPa
Type of protection against electric shock:	Internal power supply

Degree of protection against electric shock:	Device Classification Type BF applied part
Ingress Protection rating:	IP22
Expected service life:	5 years

Main Performance:

Measurement Range: 32°C – 42.9°C (89.6°F – 109.2°F)

Accuracy: $\pm 0.2^{\circ}\text{C}$ / $\pm 0.4^{\circ}\text{F}$

The clinical accuracy or clinical deviation of the medical infrared forehead thermometer has been clinically verified and can be obtained by contacting the importer of the product in your country.

X. BATTERY REPLACEMENT

1. Press the battery compartment button (9) to eject the battery compartment (see the illustration).
2. Remove the used batteries and replace them with 2 new ones (alkaline batteries 1.5V AAA). Always place the batteries correctly according to their polarity (+/-). An error with this action can lead to damage of the device and warranty canceling. Never use rechargeable batteries. Use only single use batteries.
3. Close the battery compartment properly.

IMPORTANT BATTERY GUIDELINES:

1. An adult must handle the battery.
2. Do not allow children to play with the batteries. Risk of swallowing! Small children can swallow the batteries and to suffocate. Store the batteries at place where children cannot reach!
3. Always remove exhausted batteries and do not throw them away together with the household waste, but in the intended places. They are recyclable.
4. Use batteries only AAA size.
5. Alkaline batteries are recommended.
6. When you change with new batteries, always change all batteries. Children must not be present when changing the batteries.
7. Always place the batteries correctly according to their polarity (+/-), marked on the bottom of the battery compartment. If they are not correctly placed, this can lead to short circuit or electric shock.
8. Do not mix old batteries with new ones.
9. Do not mix alkaline, standard (Carbon-Zinc) or rechargeable (Nickel- Cadmium) batteries.

10. Do not put the batteries in fire, because they can explode or leak.
11. Keep the batteries from direct heating sources.
12. If you will not use the device for a long time, take the batteries out of the battery compartment. This way you will prevent the device from battery leaking or corrosion.
13. If when changing the batteries, you find that the batteries leak, put safety gloves on and clean the battery compartment with dry towel.
14. Non-rechargeable batteries must not be charged.
15. Do not disassemble, open or shatter the batteries.
16. If your skin or eyes come into contact with battery fluid, rinse the affected area with water and seek medical assistance.
17. If the batteries have leaked, use gloves to remove the batteries and clean the battery compartment with dry cloth.
18. Protect batteries from excessive heat.

XI. CLEANING, DESINFECTION AND MAINTENANCE INSTRUCTION

CLEANING AND DESINFECTION:

CLEANING OF THE PRODUCT

1. Before cleaning the device, take out the batteries. Wipe the body of the thermometer with clean soft cloth and clean the front glass with dry cotton swab.
2. You can clean the thermometer with clean soft wet cloth and soft neutral cleaning agent. Wipe dry immediately afterwards.
3. Do not wash the thermometer under running water. Do not immerse it in water. Risk of electric shock.
4. Do not use aggressive cleaning detergents, as bleach, those with abrasive particles or on alcohol basis.
5. It is recommended to clean the thermometer once a week. Cleaning should not be longer than 3 minutes. You should not clean the product more than 2 times a day. After disinfection wait at least 10 minutes before the next use.
6. After cleaning with medical alcohol, you need to wait for 5 minutes before next – time use. Otherwise, the measurement results may be incorrect.

DESINFECTION OF THE PRODUCT

1. Wipe the body of the thermometer with a clean, soft cloth dipped with little disinfectant (70% isopropyl alcohol solution or 3% sodium hypochlorite sol) or 75% medical alcohol, then wipe it dry immediately.

MAINTENANCE

- Do not expose the thermometer on sunlight or water. Their influence will damage the product.
- If you will not use the device for long time, take the batteries out, to avoid leaking

or corrosion.

- Always store the thermometer with the cap (6).
- Do not store the device in high or low temperatures or high humidity (look technical characteristics), in direct sunlight, close to direct heating sources or dusty places. Otherwise, measurement deviations can be caused. Store the product in a dry, well-ventilated place away from direct sunlight. Make sure that the storage conditions meet the requirements.
- Do not store in places where the thermometer will be tilted or exposed to shock and vibration.
- To ensure safe use check the product regularly for damages – broken lenses, broken parts of the body, if the probe is dirty or not. Stop using the device if you find any damage.
- If you will not use the device for a longer period of time, you should clean it at least once a year.
- After each use clean the product as described in CLEANING AND DESINFECTION.

XII. TROUBLESHOOTING

Display message Problem	Possible reason	Resolvent
The thermom- eter cannot be switched on.	The battery has no charge.	Replace with new batteries.
	Battery polarity is not cor- rect.	Place the batteries according to the correct polarity as shown on the battery cover.
	The thermometer is dam- aged.	Contact the shop from which you purchased the product.
“Er1” is displayed when powered on	The ambient temperature is higher than 40°C (104.00°F) or lower than 10°C (50.0°F).	Operate the thermometer only between the specified temper- ature ranges 10°C (50.0°F) ~ 40.0°C (104.00°F).
The measured temperature is too low com- pared to normal temperature range.	Measurement area is away from the thermometer. The distance of the probe and the measuring target is too far.	Measure according to instruc- tions. Measure the temperature on recommended distance (around 0-3 cm, 0-1.18 inches).
	Temperature measuring probe (the lens of the probe) or temperature measuring area is dirty.	Measure after removing dirt.
The measured temperature is too high compared to normal tempera- ture range	The probe is damaged.	Contact the shop from which you purchased the product. Do not try to repair the product by themselves. Otherwise, the guarantee will be void.

The measured temperature is too high compared to normal temperature range	The thermometer is moved from a cold environment to the measuring environment and has not stayed for more than 30 minutes.	Before measuring temperature, the user and the thermometer should stay in the same room (premises) with constant ambient conditions for at least 30 minutes.
---	--	--

XIII. WARRANTY AND AFTER-WARRANTY SERVICE

The digital thermometer is made for intensive use; its durability is guaranteed for 1000 measurements. Any damage resulting from improper use, the batteries and the packaging are not covered by the warranty.

XIV. COMPLIANCE WITH EUROPEAN LEGISLATION AND STANDARDS

WARNING!

The infrared thermometer needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the **ACCOMPANYING DOCUMENTS**.

WARNING!

- Do not bring the thermometer close to active high-frequency surgical equipment and RF-shielded MRI rooms where the intensity of electromagnetic interference is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided as it may result in malfunction. If such use is necessary, this equipment and the other equipment should be observed to verify normal operation.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and may result in malfunction.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Failure to do so may result in deterioration of the performance of this equipment.

Note

If users or patients have experienced any serious incident related to the device, contact the retailer from whom you purchased the product or the importer.

A. APPENDIX

Guidance and manufacturer's declaration – electromagnetic emissions – for ALL EQUIPMENT AND SYSTEMS.

Guidance and manufacturer's declaration - electromagnetic emissions	
The Infrared Thermometer AOJ-20F is intended for use in the electromagnetic environment specified below. The customer or the user of the AOJ-20F should ensure that it is used in such an environment.	
Emissions test	Compliance
RF emissions CISPR 11	GROUP 1
RF emissions CISPR 11	CLASS B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable

B. ELECTROMAGNETIC IMMUNITY

Guidance and manufacturer's declaration - electromagnetic emissions		
The Infrared Thermometer AOJ-20F is intended for use in the electromagnetic environment specified below. The customer or the user of AOJ-20F should ensure that it is used in such an environment.		
Immunity test	IEC 60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/ burst IEC 61000-4-4	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable
Power frequency of the magnetic field IEC 61000-4-8	50Hz/ 60Hz	50Hz/ 60Hz
Conducted RF IEC61000-4-6	Not applicable	Not applicable
Radiated RF IEC61000-4-3	10 V/ m 80 MHz – 2,7 GHz 80 % AM at 2Hz	10 V/ m 80 MHz – 2,7 GHz 80 % AM at 2Hz
Note: UT is the a.c. mains voltage prior to application of the test level.		

C. ELECTROMAGNETIC IMMUNITY

Guidance and manufacturer's declaration – electromagnetic immunity				
Radiated RD IEC 61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Modulation	IEC 6061-1-2 Test Level (A/m)	Compliance level (A/m)
	30 kHz	CW	8	8
	134.2 kHz	Pulse modulation 2.1 kHz	65	65
	13.56 MHz	Pulse modulation 50 kHz	7.5	7.5

XIV. SYMBOLS

Symbol	Explanation
	Trade mark
	Medical device
	The CE labelling certifies that the product complies with the essential requirements of EU regulations on medical device.
	Manufacturer
	Authorized Representative in the European Community
	Importer
	Made in country
	Item Number
	Warning
	IEC 60417-5333, Type BF applied part
	IEC 60417-5031, Direct current
IP22	2 Protected against solid foreign objects of 12,5 mm Ø and greater; 2 If you keep the thermometer in 15 degree angle, it still can prevent the water drop.
	Directing to the Instruction/ Manual
	DISPOSING: Do not throw this product with household waste. Collect separately, because special disposal of those kind of waste is required.
	Unique device identifier
	Lot number
	Specific serial number
	Date of Manufacture
	Use-by date

RoHS	RoHS symbol
	Keep dry
	Fragile

Origin China

Shenzhen AOJ Medical Technology Co., Ltd.

Add: Room 301&4F, Block A, Building A, Jingfa Intelligent Manufacturing Park, Xiaweyuan, Gushu Community, Xixiang Street, Bao'an District, 518126 Shenzhen, CHINA

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EN WARRANTY: The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions. The warranty shall not apply in the case of damage caused by improper use, wear or accidental events. For the duration of the warranty on conformity defects please refer to the specific provisions of applicable national standards in the country of purchase, where provided.

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Share Info GmbH
Moni Trade LTD
Manufacturer and Importer: Moni Trade Ltd.