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Shenzhen AOJ Medical Technology Co., Ltd.

Add: Room 301&4F, Block A, Building A, Jingfa Intelligent Manufacturing Park, Xiawei-yuan, Gushu Community, Xixiang Street, Bao'an District, 518126 Shenzhen, CHINA

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THE MANUFACTURER KEEPS THEIR RIGHTS TO CHANGE SPECIFICATIONS OF THE PRODUCT WITHOUT PRE-WARNING!

IMPORTANT! READ CAREFULLY AND KEEP FOR FUTURE REFERENCE! USE THIS DEVICE EXACTLY AS DESCRIBED IN THE MANUAL AND ONLY IF IT IS IN PROPER WORKING CONDITION.

I. INTENDED USE

The appliance is infrared thermometer is intended for measuring the body temperature of both children and adults. The temperature can be measured either by pointing the thermometer at the child's or adult's forehead or by auricular measurement (outer ear canal). Measuring the temperature by the auricular method is suitable for children above 3 months of age. The device can be used in home care environment. The product is tested according to and complies with the requirements of ISO 81060-2-56, ISO10993-5, ISO10993-10 and ISO10993-23.

II. INCLUDED IN THE PACKAGE

Before use, carefully open the package and check that all parts are present as listed below and have not been damaged during shipment. If there is any damage or problem related to operation, contact the dealer from whom you purchased the product.

INCLUDED IN THE PACKAGE:

Infrared thermometer – 1 pc;

Pouch – 1 pc;

User manual – 1 pc;

Batteries AAA 1.5V – 2 pcs.

III. RULES AND RECOMMENDATIONS FOR SAFE USAGE

1. Only use this thermometer once you have read and understood the instruction for use.
2. Only use this thermometer for the specific purpose described in the instruction manual. Any misuse will void the warranty.
3. Do not place the thermometer near the reach of children as this may lead to damaging the product.
4. Children must not play with the device.
5. Device can be used for personal needs.
6. Do not use the thermometer to measure the temperature of newborn babies or for continuous temperature monitoring purposes.
7. This thermometer is not intended for preterm babies. This unit is not intended to interpret hypothermic temperatures.
8. Do not measure the temperature of the baby during or right after nursing.
9. The patients must not drink, eat or to be physically active before or during using the thermometer.
10. Using this thermometer is no substitute for consulting a doctor. Do not self-diagnose or self-medicate based on the measurements without consulting your doctor. In particular, do not start taking any new medication or change the type and/ or dosage of any existing medication without prior approval. Consult your doctor if you have symptoms such as vomiting, diarrhea, dehydration, loss of appetite and fatigue, seizures, shivering, stiff neck, pain when urinating, etc., even in the absence

of fever.

11. This thermometer must be operated by adults and be placed out of the reach of infants.

12. Warning! Small parts. Risk of swallowing or choking. Make sure that children do not use the device unsupervised.

13. If this thermometer is stored in a cooler or warmer location than where it is being used, put it in the operation place for at least 30 minutes before taking the measurement.

14. After physical activity, a bath, or spending time outdoors, wait approximately 30 minutes before you measure the temperature on your forehead.

15. To ensure an accurate measurement, you should only use this thermometer indoors when the ambient temperature range is between 10°C (50°F) and 40°C (104°F). Do not use the thermometer in environment with high humidity (>95% RH).

16. When the temperature in the room is below 10°C, place the device at a place with temperature between 10°C and 40°C for at least 1 hour. When the temperature in the room is above 40°C, place the device at a place with temperature between 10°C and 40°C for at least 2 hours. Extreme temperatures, humidity and altitude may lead to inaccurate measurements.

17. Before and after each use, the thermometer should be cleaned according to the requirement of this manual to avoid possible bacterial infections.

18. Do not forcibly insert the probe into the ear canal, otherwise you may injure the ear canal. If the ear canal is dirty, first clean it with cotton swab and then measure the temperature.

19. Do not measure your temperature auricularly if you have an ear infection - otitis.

20. Do not expose this thermometer to electric shocks.

21. Do not expose the thermometer on extreme temperature conditions, bigger than 40°C and less than 10°C.

22. Do not use the device in relative humidity, higher than 95%.

23. Do not put this product in a place with direct sunlight, high temperature, high humidity, excessive dust, or the place where it may contact corrosive gas or liquid. Do not use this product in the above environment.

24. Never drop the device.

25. DO NOT TOUCH THE GLASS OF THE INFRARED SENSOR WITH YOUR FINGERS AS IT IS EXTREMELY FRAGILE AND COULD CAUSE INJURY.

26. Do not use this thermometer if it is damaged.

27. Do not open up this thermometer. It contains small parts which might be swallowed by children.

28. Never immerse the device in water or other cleaning liquids, and do not let any liquid enters the interior of this product.

29. If you will not use the device for a time longer than 2 months, take the batteries out from the battery compartment. This way you will prevent the device from battery leaking or corrosion.

30. Keep the thermometer away from electromagnetic fields produced by objects such as microwave ovens, arc welders or induction heaters.

31. Portable and mobile radio frequency communication devices may affect the operation of the infrared thermometer.

32. The device should no longer be reused when it reaches the expected service life (5 years, from the date of production) of the device. Electrical components and batteries have heavy metal or chemicals which may do harm to environment, follow the local law and regulation on its disposal.

33. If a problem with the device appears, contact with the shop, from which you have bought it. Do not try to repair the device by yourself. Otherwise, the warranty will be void:

- If an error occurs or if the thermometer is damaged. DO NOT attempt to repair the thermometer yourself, as this will void the warranty.
- Only have repairs carried out by authorized service partners.
- If you need assistance in setting up, using, or maintaining the thermometer.
- If you have unexpected operation or events to report.

34. The replaceable parts of the thermometer are the 2 batteries AAA 1.5V. Replace them as soon as they have no charge. Do not use batteries with different voltages or different specifications. The product does not have any other replaceable parts and materials.

35. Even in the absence of fever, those who exhibit a normal temperature may still need medical attention. People who are on antibiotics, analgesics or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.

36. A rise in temperature may be a sign of a serious illness, especially in elderly individuals, those who are frail, have weakened immune systems, or in newborns and infants. Seek professional advice immediately when there is an increase in temperature and if you are measuring the temperature of:

- People over 60 years of age (Fever may be reduced or even absent in elderly patients).
- People with diabetes or weakened immune systems (e.g. HIV positive, cancer patients, undergoing chemotherapy, on long-term steroid treatment, or with a splenectomy).
- Bedridden patients (e.g. those in nursing homes, with heart attacks, or chronic illnesses).
- Transplant patient (e.g. kidney, heart, lung, or liver transplants).

IV. WORKING PRINCIPLE

All objects, solid, liquid or gas emit energy by radiation. The intensity of this energy depends of the temperature in the object. The infrared thermometer Easy Check is therefore able to measure the human's temperature, only by the energy that person emits. This measurement is possible, thanks to an external temperature probe on the device, which permanently analyses and register the temperature of the environment. Therefore, immediately after the thermometer is closed to the body and the radiation sensor is activated, the temperature is immediately measured through detecting infrared heat, caused by arterial blood flow. The temperature of the body can be measured also without any interference from the heating of the surrounding environment. The temperature can be measured and by placing the sensor in the ear. The product is with alarm for high temperature, front light, mute button, function to change the measurement unit °C and °F, automatic turning off to save energy.

V. DEVICE DESCRIPTION


FIGURE A: Device

1. LCD display screen.
2. Mode button/ Memory (switch for body temperature/ room temperature/ object).
3. Measure/ Mute-unmute/ Turn on button.
4. Infrared sensor/ probe (measuring the body temperature in the ear must be made without the cover).
5. Battery compartment.
6. Probe cover (measuring the body temperature on the forehead is made with the probe cover on).

FIGURE B: Display

1. Adult forehead temperature mode (above 12 years old).
2. Child forehead temperature mode (below 12 years old).
3. Ear temperature mode (measuring the temperature in the ear of an adult and children above 3 months of age).
4. Mode for measuring the temperature of an object/ room temperature/ adjust mode.
5. Celsius degree °C. Fahrenheit degree °F.
6. Previous temperature value.
7. Low battery indicator.
8. Mute/Unmute indicator.
9. Low battery indicator.
10. Current temperature value.

VI. HOW TO USE THE THERMOMETER:

1. To avoid incorrect measurements:


- Make sure that the device is used only indoors and there is no strong draft.
- Make sure that the patient has not experienced strong emotions or performed intense physical activity before taking the temperature.
- If the thermometer was stored or used in a room with lower or higher temperature than the temperature of the room you intend to use it, you should leave the device for 30 minutes in the new room before use.
- If the patient has been in an environment with a lower or higher temperature than that of the room in which their temperature will be measured, they must

wait at least 10 minutes before measuring.

- Do not hold the unit for too long as it is sensitive to heat.
- The device has passed the necessary tests and is safe to use, and provides accurate results when used correctly according to the recommendations in the instructions.

2. Measuring the forehead temperature:

- Place the batteries in the compartment before first time use.
- Press the button for Measuring/ Mute - unmute (3), to power on the thermometer. Do not hold the button for more than 1 second. Press the Mode button (2) and choose from  and . Point the probe of the thermometer to the center of the forehead at a distance of 0 – 3 cm (0-1.18 inch).
- Press the button Measure/ Mute -unmute for not more than 1 second. You will hear indicating sound and the temperature will be displayed on the screen.

NOTE: Measuring the temperature on the forehead is an indicative reading. The measured forehead temperature may vary with 1 °F / 0.5 °C from your actual body temperature. Keep in mind the factors that can influence the accuracy as described in ADVICES WHEN MEASURING THE TEMPERATURE AND RULES AND RECOMMENDATIONS FOR SAFE USE.

If the forehead is covered with hair, sweat, or dirt, clean it to ensure an accurate temperature measurement.

Always check if the lens is clean.

Ensure that the thermometer has been in the same room as the patient for at least 30 minutes before taking the measurement.

3. Auricular temperature measurement:

- Press the button Measure/ Mute-unmute (3) on power on the thermometer.
- Press the Mode (2) button and choose mode  and remove the probe cover (if the cover is magnetic, the measurement from forehead to ear is switched automatically) and place the probe carefully in the ear canal.
- Press the measurement/ Mute-unmute (3) button for not more than 1 second. You will hear an indicating sound and the temperature will be displayed on the screen.

• Press the measurement/ Mute-unmute (3) button for not more than 1 second. You will hear an indicating sound and the temperature will be displayed on the screen.

Make sure that the ear is clean (no earwax or any inflammation).

The result from measuring the temperature in the right ear can be different from the one, when measuring the temperature in the left ear. We recommend measuring the temperature always in the same ear.

FOR CHILDREN UNDER ONE YEAR OLD: As the ear canal is slightly curved, you have to pull the ear slightly backwards before inserting the sensor tip. This is important so that the sensor tip can be pointed directly at the ear drum.

FOR CHILDREN ABOVE ONE YEAR OLD AND ADULTS: As the ear canal is slightly curved, you have to pull the ear slightly up and backwards before inserting the sensor tip. This is important so that the sensor tip can be pointed directly at the ear drum.

WARNING! Do not forcibly insert the probe into the ear canal; otherwise, you may injure the ear canal.

WARNING! Be careful when measuring the ear temperature of a child, as a child's earhole is small.

4. After measuring the temperature:

- After the temperature is measured remove the thermometer from the forehead/ ear and observe the result.
- After each measurement you can enter the memory mode and check the previous measurements.
- Do not hold the thermometer for too long and is it sensitive to ambient temperature.

After each measurement clean the sensor with cotton swab and store the thermometer in a dry and well-ventilated place.

- You should wait at least 10 seconds between each measurement.
- Self-diagnosis or self-treatment is not recommended without consulting your doctor or relying solely on the measured temperature.

For proper treatment, always consult your doctor.

5. Reading the measured temperature – T stands for measured temperature:

- **In forehead or ear mode:** If the measured temperature is in the range of $32^{\circ}\text{C} \leq T \leq 37.5^{\circ}\text{C}$ ($89.6^{\circ}\text{F} \leq T \leq 99.5^{\circ}\text{F}$), the display will light up green and a single long beep will sound. If the measured temperature is in the range of $37.6^{\circ}\text{C} \leq T \leq 38.30^{\circ}\text{C}$ ($99.7^{\circ}\text{F} \leq T \leq 100.94^{\circ}\text{F}$), the display will turn orange and short beeps will sound, indicating a mild fever. If the measure temperature is in the range of $38.1^{\circ}\text{C} \leq T \leq 42.9^{\circ}\text{C}$ ($100.6^{\circ}\text{F} \leq T \leq 109.2^{\circ}\text{F}$), the display will light up red and six short beeps will sound, indicating high fever.

- **In adjust mode:** If the measured temperature is in range $0^{\circ}\text{C} \leq T \leq 100^{\circ}\text{C}$ ($32.0^{\circ}\text{F} \leq T \leq 199.0^{\circ}\text{F}$), the display will light up white and one long beep will sound. This mode is not for measuring the human temperature, as it is only for maintenance.

6. Switching to silent mode: When the thermometer is switched on, press and hold the Measurement/ Mute-unmute button for approximately 5 seconds to switch from unmute to mute or vice versa.

7. Changing the measurement unit – $^{\circ}\text{C}$ / $^{\circ}\text{F}$: Open the battery compartment cover and press the button for switching between $^{\circ}\text{C}$ or $^{\circ}\text{F}$ (Figure I – button B):

8. Automatic shutdown: The device will switch off automatically if you do not use it for 60 ± 10 seconds.

Warning! All memory records and settings will be lost when you replace the batteries or remove and reinsert them.

9. Changing the batteries: Place two batteries AAA correctly according to the polarity in the battery compartment.

Remove the batteries from the thermometer if you will not use it for more than 2 months.

10. IMPORTANT BATTERY GUIDELINES:

1. An adult must handle the battery.

2. Do not allow children to play with the batteries. Risk of swallowing! Small children can swallow the batteries and to suffocate. Store the batteries at place where children cannot reach!

3. Always remove exhausted batteries and do not throw them away together with the household waste, but in the intended places. They are recyclable.

4. Use batteries only AAA size.
5. Alkaline batteries are recommended.
6. When you change with new batteries, always change all batteries. Children must not be present when changing the batteries.
7. Always place the batteries correctly according to their polarity (+/-), marked on the bottom of the battery compartment. If they are not correctly placed, this can lead to short circuit or electric shock.
8. Do not mix old batteries with new ones.
9. Do not mix alkaline, standard (Carbon-Zinc) or rechargeable (Nickel- Cadmium) batteries.
10. Do not put the batteries in fire, because they can explode or leak.
11. Keep the batteries from direct heating sources.
12. If you will not use the device for a long time, take the batteries out of the battery compartment. This way you will prevent the device from battery leaking or corrosion.
13. If when changing the batteries, you find that the batteries leak, put safety gloves on and clean the battery compartment with dry towel.
14. Non-rechargeable batteries must not be charged.
15. Do not disassemble, open or shatter the batteries.
16. If your skin or eyes come into contact with battery fluid, rinse the affected area with water and seek medical assistance.
17. If the batteries have leaked, use gloves to remove the batteries and clean the battery compartment with dry cloth.
18. Protect batteries from excessive heat.

VII. ADVICE WHEN MEASURING THE TEMPERATURE:

1. Always measure the temperature at the same place in order to receive accurate results.
2. A child's body temperature can vary. Note that this device measures approximately 0.5°C (0.9°F) lower than a rectal digital thermometer.
3. Certain external factors may affect the accurate measurement of ear temperature:
 - If you have been lying on one ear or if the ears were covered.
 - If you have been exposed to very high or very low temperatures.
 - If you have recently been swimming or bathing.
 - In such cases, wait at least 20 minutes before measuring the temperature. If one ear is being treated with drops or other medications, measure the temperature in the other ear.
4. Holding the thermometer in your hand for an extended period may cause it to warm up and result in inaccurate readings.
5. Both the patient and the thermometer should remain at a constant room temperature for at least 30 minutes before measurement.
6. Clean the forehead of any dirt, hair, or sweat before placing the thermometer sensor on it. Wait 10 minutes after cleaning before measuring the temperature.
7. Use an alcohol-soaked cotton pad to gently clean the sensor, and wait 5 minutes before measuring another patient's temperature. Wiping the forehead with a warm or cold cloth can affect the measurement result. We recommend waiting at least 10 minutes before taking the temperature.

8. In the following cases, it is advisable to take 3 to 5 measurements at the same spot and use the highest value as the correct result:

- In newborns during the first 100 days.
- Children under three years of age with weakened immune systems, for whom the presence or absence of a fever is critically important.
- When the user is using the thermometer for the first time, until they become familiar with the device and begin to obtain consistent results.

VIII. CLEANING, DESINFECTION AND MAINTENANCE INSTRUCTION

CLEANING AND DESINFECTION: CLEANING OF THE PRODUCT

1. Before cleaning the device, take out the batteries. Wipe the body of the thermometer with clean soft cloth and clean the front glass with dry cotton swab.
2. You can clean the thermometer with clean soft wet cloth and soft neutral cleaning agent. Wipe dry immediately afterwards.
3. Do not wash the thermometer under running water. Do not immerse it in water. Risk of electric shock.
4. Do not clean the prob lens with a paper towel, as this may scratch the glass and lead to incorrect results.
5. Do not use aggressive cleaning detergents, as bleach, those with abrasive particles or on alcohol basis.
6. It is recommended to clean the thermometer once a week. Cleaning should not be longer than 3 minutes. You should not clean the product more than 2 times a day. After disinfection wait at least 10 minutes before the next use.
7. After cleaning with medical alcohol, you need to wait for 5 minutes before next – time use. Otherwise, the measurement results may be incorrect.

8. Be careful when cleaning not to scratch the LCD screen.
9. Cleaning and disinfection must be done in the following conditions: +10°C ~ 40°C with relative humidity 15-95% RH (without condensing), atmospheric pressure 10kPa ~ 106kPa.

DESINFECTION OF THE PRODUCT

1. Wipe the body of the thermometer with a clean, soft cloth dipped with little disinfectant (70% isopropyl alcohol solution or 3% sodium hypochlorite sol) or 75% medical alcohol, then wipe it dry immediately.

MAINTENANCE

- Do not expose the thermometer on sunlight or water. Their influence will damage the product.
- If you will not use the device for long time, take the batteries out, to avoid leaking or corrosion.
- Always store the thermometer with the cap (6).
- Do not store the device in high or low temperatures or high humidity (look technical characteristics), in direct sunlight, close to direct heating sources or dusty places. Otherwise, measurement deviations can be caused. Store the product in a dry, well-ventilated place away from direct sunlight. Make sure that the storage conditions meet the requirements.
- Do not store in places where the thermometer will be tilted or exposed to shock and vibration.

- To ensure safe use check the product regularly for damages – broken lenses, broken parts of the body, if the probe is dirty or not. Stop using the device if you find any damage.
- If you will not use the device for a longer period of time, you should clean it at least once a year.
- After each use clean the product as described in CLEANING AND DESINFECTION.
- If the device is damaged and cannot be repaired, dispose of it at designated collection points in accordance with local environmental protection regulations.

IX. PRODUCT TECHNICAL SPECIFICATIONS

Model:	AOJ-20A
Battery-powered mode/ Batteries	DC 3V/ 2 batteries 1.5 V AAA
Display mode:	LCD (4 color backlight display-white, green,red,orange)
Display temperature resolution:	0.1°C/0.1°F
Working mode:	Continuous operation
Measurement time:	Short recall time - around 1 second
Measuring range:	Forehead/ Ear temperature: 32°C – 42.9°C (89.6°F – 109.2°F) ±0,2°C /±0,4°F.
	Adjusted temperature (room/object temperature): 0.0°C-100.0°C (32.0°F-199.0°F) ±1.0°C/±2.0°F.
Display resolution:	0.1°C
Automatic shutdown:	The thermometer will shut down automatically after 60±10 seconds if is not being used.
Memory data:	40 measurements
Product size:	14.0 cm x 3.5 cm x 4.0 cm
Weight (g):	about 60 g (without batteries)
Normal operating conditions (Using environment):	- Environment temperature: from +10°C to +40°C (+50°F - +104°F) - Relative humidity: 15%~95% RH - Atmospheric pressure: 70kPa~106kPa
Storage and transportation condition (storage environment):	- Environment temperature: -20°C~+55°C (-4°F-131°F) - Relative humidity: 0~95% RH - Atmospheric pressure: 50kPa~106kPa
Type of protection against electric shock:	Internal power supply

Degree of protection against electric shock:	Device Classification Type BF applied part
Ingress Protection rating:	IP22
Expected service life:	5 years

Main Performance:

Measurement Range: 32°C – 42.9°C (89.6°F – 109.2°F)

Accuracy: $\pm 0.2^{\circ}\text{C}$ / $\pm 0.4^{\circ}\text{F}$

The clinical accuracy or clinical deviation of the medical infrared forehead thermometer has been clinically verified and can be obtained by contacting the importer of the product in your country.

The infrared thermometer uses Bluetooth 4.2 technology with transmission and reception frequencies ranging from 2402 to 2480 MHz, GFSK modulation type, and an effective radiated power of 2.79 dBm.

X. TROUBLESHOOTING

Display message Problem	Possible reason	Resolvent
The thermometer cannot be switched on.	The battery has no charge.	Replace with new batteries.
	Battery polarity is not correct.	Place the batteries according to the correct polarity as shown on the battery cover.
	The thermometer is damaged.	Contact the shop from which you purchased the product.
The measured temperature is too low compared to normal temperature range.	Measurement area is away from the thermometer. The distance of the probe and the measuring target is too far.	Measure according to instructions. Measure the temperature on recommended distance (around 0-3 cm, 0-1.18 inches).
	Temperature measuring probe (the lens of the probe) or temperature measuring area is dirty.	Clean thoroughly and repeat the measurement.

Display message Problem	Possible reason	Resolvent
The measured temperature is too low compared to normal temperature range.	The patient was in an environment with too low temperature.	The patient must wait at least 30 minutes in a room with warmer temperature before measuring their temperature.
The measured temperature is too high compared to normal temperature range	The probe is damaged.	Contact the shop from which you purchased the product. Do not try to repair the product by themselves. Otherwise, the guarantee will be void.
	The thermometer is moved from a cold environment to the measuring environment and has not stayed for more than 30 minutes.	Before measuring temperature, the user and the thermometer should stay in the same room (premises) with constant ambient conditions for at least 30 minutes.
	The patient was in an environment with too high temperature.	The patient must wait at least 30 minutes in a room with cooler temperature before measuring their temperature.
Er1" is displayed when powered on	The ambient temperature is higher than 40°C (104.00°F) or lower than 10°C (50.0°F).	Operate the thermometer only between the specified temperature ranges 10°C (50.0°F) ~ 40.0°C (104.00°F).
"ErC" is displayed when powered on	Memory error or device calibration is incomplete. The device will automatically turn off after 3 seconds.	If the display lights up red and a long beep sounds, contact the retailer where you purchased the product.

DISPLAY OPERATIONS

SCREEN DISPLAY	DESCRIPTION	LCD INDICATOR STATUS
DISPLAY SHOWS THE MESSAGE "HI"		
	In this case the temperature is either over the chosen range or over 42,9 °C (109,2 °F) in „Ear temperature“ or „Forehead temperature“ mode. In adjust mode the measured temperature is >100°C (199°F).	One long beep and LCD indicator in red
DISPLAY SHOWS MESSAGE "LO"		
	In this case, recorded temperature is either under chosen range or under 32 °C (89,6 °F) in „Ear temperature“ or „Forehead temperature“ mode. In adjust mode the measured temperature is >0°C (32°F).	One long beep and LCD indicator in red
BATTERY LEVEL INDICATOR		
	The battery charge is between 2.4V – 2.7V.	The battery icon will flash, indicating that the battery is low. However, you will still be able to measure the temperature.
	The battery charge is lower than 2.4V.	The battery icon will flash, indicating that the battery is low. The device will automatically turn off after 3 seconds. Replace the depleted batteries with new ones.

XI. WARRANTY AND AFTER-WARRANTY SERVICE

The digital thermometer is made for intensive use; its durability is guaranteed for 1000 measurements. Any damage resulting from improper use, the batteries and the packaging are not covered by the warranty.

XII. COMPLIANCE WITH EUROPEAN LEGISLATION AND STANDARDS

WARNING!

The infrared thermometer needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the **ACCOMPANYING DOCUMENTS**.

WARNING!

- Do not bring the thermometer close to active high-frequency surgical equipment and RF-shielded MRI rooms where the intensity of electromagnetic interference is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided as it may result in malfunction. If such use is necessary, this equipment and the other equipment should be observed to verify normal operation.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and may result in malfunction.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Failure to do so may result in deterioration of the performance of this equipment.

Note

If users or patients have experienced any serious incident related to the device, contact the retailer from whom you purchased the product or the importer.

A. APPENDIX

Guidance and manufacturer's declaration – electromagnetic emissions – for ALL EQUIPMENT AND SYSTEMS.

Guidance and manufacturer's declaration - electromagnetic emissions	
The Infrared Thermometer AOJ-20A is intended for use in the electromagnetic environment specified below. The customer or the user of the AOJ-20A should ensure that it is used in such an environment.	
Emissions test	Compliance
RF emissions CISPR 11	GROUP 1
RF emissions CISPR 11	CLASS B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable

B. ELECTROMAGNETIC IMMUNITY

Guidance and manufacturer's declaration - electromagnetic emissions		
The Infrared Thermometer AOJ-20A is intended for use in the electromagnetic environment specified below. The customer or the user of AOJ-20A should ensure that it is used in such an environment.		
Immunity test	IEC 60601-1-2 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/ burst IEC 61000-4-4	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable
Power frequency of the magnetic field IEC 61000-4-8	50Hz/ 60Hz	50Hz/ 60Hz
Conducted RF IEC61000-4-6	Not applicable	Not applicable
Radiated RF IEC61000-4-3	10 V/ m 80 MHz – 2,7 GHz 80 % AM at 2Hz	10 V/ m 80 MHz – 2,7 GHz 80 % AM at 2Hz
Note: UT is the a.c. mains voltage prior to application of the test level.		

C. ELECTROMAGNETIC IMMUNITY

Guidance and manufacturer's declaration – electromagnetic immunity				
Radiated RD IEC 61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Modulation	IEC 60601-1-2 Test Level (A/m)	Compliance level (A/m)
	30 kHz	CW	8	8
	134.2 kHz	Pulse modulation 2.1 kHz	65	65
	13.56 MHz	Pulse modulation 50 kHz	7.5	7.5

C. ELECTROMAGNETIC IMMUNITY

[illegible]

XIII. SYMBOLS

Symbol	Explanation
	Trade mark
	Medical device
	The CE labelling certifies that the product complies with the essential requirements of EU regulations on medical device.
	Manufacturer
	Authorized Representative in the European Community
	Importer
	Made in country
	Item Number
	Warning
	IEC 60417-5333, Type BF applied part
	IEC 60417-5031, Direct current
IP22	2 Protected against solid foreign objects of 12,5 mm Ø and greater; 2 If you keep the thermometer in 15 degree angle, it still can prevent the water drop.
	Directing to the Instruction/ Manual
	DISPOSING: Do not throw this product with household waste. Collect separately, because special disposal of those kind of waste is required.
	Unique device identifier
	Lot number
	Specific serial number
	Date of Manufacture
	Use-by date

RoHS	RoHS symbol
	Keep dry
	Fragile

Origin China

Shenzhen AOJ Medical Technology Co., Ltd.

Add: Room 301&4F, Block A, Building A, Jingfa Intelligent Manufacturing Park, Xiaweyuan, Gushu Community, Xixiang Street, Bao'an District, 518126 Shenzhen, CHINA

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EN WARRANTY: The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions. The warranty shall not apply in the case of damage caused by improper use, wear or accidental events. For the duration of the warranty on conformity defects please refer to the specific provisions of applicable national standards in the country of purchase, where provided.

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Share Info GmbH
Moni Trade LTD
Manufacturer and Importer: Moni Trade Ltd.